

# Max 2M™

MRI contrast media injector



Contrast-enhanced  
mammography



**Ulrich**  
medical

*dedicated to you.*

# Key new functions

*Your workflow has never been easier*



## Contrast-enhanced mammography

### Fast operation

Patient-side injection start directly on the injector.

### Space-saving

Compact device base ensures sufficient space in the mammography room.

### Reproducibility

Using a contrast media injector provides a reproducible and continuous flow of contrast media

### Improved image quality

CM bolus is administered accurately, which then increases the quality of the imaging.

### Flexibility

Cordless device concept for free and flexible positioning of the injector in the mammography room.

## Keep-Vein-Open

- Quick and easy access to the KVO protocol
- The KVO protocol can be selected via the program list (fixed protocol) on the touch terminal
- KVO can be selected directly via the main menu on the injector display
- “Pause time” can be set from 15-75 s; injection: 0.3 ml with a flow of 0.3 ml/s per set pause time

**NEW**


Expanded  
KVO functions






## Hardware Features

*User-friendly.  
Simple. Efficient.*



### **Set-up made easy**

- Insert Easy-Click-Cassette (Click-Click)
- Connect patient tubing
- Attach media containers
- Done!



### **Patient tubing changeover in seconds**


- Detach patient tubing
- Connect new patient tubing
- Done!



### **Powerful touch terminal for convenient working**

- Wireless connection to injector
- Fast and uninterrupted connection via WLAN

# Innovative technology made in Germany



## Patient tubing

- Patient-side Luer lock connector
- Two check valves

## Media connection points

- Two media connection points for standard media containers
- Status display

## Injector display

- Intuitive and simple user interface
- Touch display directly on the injector

## Easy-Click-Cassette

- 24 h use
- Air and particle filters
- Connector with touch protection (SafeConnect)

## Cordless operation

- Cordless operation
- 24 h battery output
- Lithium ion battery





# Reliability and precision. For over 110 years.



## Simple to use

- Set-up just once daily
- Rapid patient changeover
- Optimal contrast media management



## Hygiene & patient safety

- Closed system
- Direct injection from the original media containers
- SafeConnect with touch protection
- Unique, automatic, patient-side air detectors protect against air injection
- Easy-Click-Cassette approved for 24-hour use<sup>1</sup>
- Reliable protection against retrograde contamination<sup>2</sup>



## Cost-effective and efficient

- Low cost of disposables
- Optimized contrast media management
- Ideal for high patient throughput



## Individual all-round service and support

- Extensive network of highly qualified service engineers
- Expert staff training and on-site application support
- Telephone support and assistance
- Tailored service contracts
- Fast repair and loan unit service
- Regular technical safety checks

<sup>1</sup>Source: GLP-compliant study: "Microbiological Contamination Control Study of the ulricheasy INJECT Max 3 System", Report 12046467, dated: February 28, 2018, author: UL International GmbH (Grenzenstr. 13, 88416 Ochsenhausen, Germany); study is unpublished, but available on request with nondisclosure agreement.

<sup>2</sup>Source: GLP-compliant study: "Virological evaluation of Retrograde Contamination Potential at the Pump Tubing Set (XD 8003) and the Patient Tubing Set (XD 2035) of the CT motion (XD 8000)", Report 10874251 2.1, dated: November 12, 2015, author: UL International GmbH (Grenzenstr. 13, 88416 Ochsenhausen, Germany); study is unpublished, but available on request with nondisclosure agreement. The results of the non-clinical investigation mentioned are applicable to the ulricheasyINJECT systems Max 3 and Max 2M as the physical barrier with respect to possible retrograde contamination is implemented in all of the systems by means of structurally identical two check valves.

# Technical data

<b>Dimensions (L×W×H)</b>	530 mm × 530 mm × 1,372 mm
<b>Weight</b>	Approx. 40 kg
<b>Power supply</b>	Battery
<b>Voltage supply (battery charging)</b>	Voltage rating 100-240 VAC / 50/60 Hz
<b>Media transport</b>	Roll pump
<b>Flow rate</b>	0.1 to 10 ml/s, in increments of 0.1 ml/s (MRI) 0.1 to 4 ml/s, in increments of 0.1 ml/s (contrast-enhanced mammography)
<b>Maximum injection volume</b>	400 ml/patient 1 ml – 10 ml, in increments of 0.1 ml > 10 ml – 100 ml, in increments of 0.5 ml > 100 ml – 400 ml, in increments of 1 ml
<b>Reservoir capacity</b>	Contrast media max. 1 × 200 ml (MRI), max. 1 × 250 ml (contrast-enhanced mammography) NaCl: max. 1 × 1,000 ml or max. 1 × 2,000 ml with a bag and optional media rod
<b>Adjustable pressure limit</b>	2-14 bar in increments of 0.1 bar (29-203 psi in increments of 1 psi)
<b>Air monitoring</b>	3 detectors for the media supply tubing 1 detector for the Easy-Click-Cassette and patient tubing
<b>Pressure monitoring</b>	1 pressure sensor
<b>User interfaces</b>	Terminal: 12" color TFT touch screen Injector: 4" color TFT touch screen
<b>Data import and export</b>	USB interface
<b>Data transfer from injector to terminal</b>	WLAN (802.11 b/g/n)
<b>Quantity of boli (phases) per injection program</b>	40 max.
<b>Quantity of injection programs that can be saved</b>	200 max.
<b>Pressure curve</b>	Real-time display of the injection pressure curve
<b>Software functions (standard)</b>	Display of remaining quantity NaCl backup Manual pause, time-controlled pause Timer Elapsed time Keep Vein Open (KVO)
<b>Software options</b>	Vein-check Same patient: Multiple injections per patient
<b>Injector accessories</b>	Convenience package: Media rod, waste bin and tray
<b>Disposables</b>	Patient tubing 2.5 m or 3.2 m Easy-Click-Cassette with holder s or holder l Adapter ring for small bottles

The country-specific availability of articles must be taken into consideration. ulrich medical does not sell and distribute this product in the United States.