


# cosmicMIA™

dynamic minimally invasive system



# System

- Dorsal dynamic stabilization, from the thoracic to the sacral spine
- Non-fusion technique – preservation of mobility in the segment
- Screw with hinged joint
- Cannulated screw with self-cutting thread
- Bonit®coating
- Minimally invasive instrumentation



3 screw diameters



5.5 mm

6.5 mm

7.5 mm



## Advantages

- Reliable pedicle screw technique
- Preservation of mobility in the segment
- Rotation-stable joint design
- Sterile packed implants



# Components

Implants	Product number
cosmicMIA™ screw, Ø 5.5 mm, lengths 30 mm – 55 mm	CS 2600-5-30 – CS 2600-5-55
cosmicMIA™ screw, Ø 6.5 mm, lengths 30 mm – 55 mm	CS 2600-6-30 – CS 2600-6-55
cosmicMIA™ screw, Ø 7.5 mm, lengths 30 mm – 55 mm	CS 2600-7-30 – CS 2600-7-55
cosmicMIA™ locking screw, blue	CS 2605
cosmicMIA™ rod, blue, Ø 6.0 mm, lengths 40 mm – 200 mm	CS 2608-040 – CS 2608-200

patented  
or/and  
pat. pend.

Over a Century  
of Innovation

**ulrich**  
medical

ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany  
Phone: +49 (0)731 9654-0 | Fax: +49 (0)731 9654-2705  
spine@ulrichmedical.com | www.ulrichmedical.com