

obelisc™

vertebral body replacement





Ulrich
medical

dedicated to you.

System

- Expandable vertebral body replacement for the anterior thoracic and lumbar spine
- Continuously expandable for adjustment to the individual anatomy
- Round, oval and rectangular end pieces for flexible adjustment to the patient
- Only one instrument for implantation, distraction and securing
- Can be used in a minimally invasive manner
- All implants are also available sterile packed




Different shapes and sizes of end pieces

∠ 0° | 5° | 10° | 15° | 20°



Different angled end pieces




Center piece with fixed or variable end pieces



Advantages

- Smooth, in situ distraction
- Multiple precise configuration options
- Allows all surgical approaches



Components

Implants	Height	Angle	Prod. No.
obelisc locking screw			CS 2901
obelisc center piece, Ø 20 mm	17 – 23 mm		CS 2920-17
	20 – 28 mm		CS 2920-20
	23 – 31 mm		CS 2920-23
	27 – 38 mm		CS 2920-27
	32 – 47 mm		CS 2920-32
	40 – 62 mm		CS 2920-40
	53 – 87 mm		CS 2920-53
	76 – 132 mm		CS 2920-76
obelisc end piece, grey, Ø 20 mm		0°/5°/10°/15°	CS 2925-00/-05/-10/-15
obelisc end piece, light blue, Ø 24 mm		0°/5°/10°/15°	CS 2926-00/-05/-10/-15
obelisc end piece, yellow, Ø 26 mm		0°/5°/10°/15°	CS 2927-00/-05/-10/-15
obelisc end piece, green, Ø 29 mm		0°/5°/10°/15°	CS 2928-00/-05/-10/-15
obelisc end piece, blue, Ø 32 mm		0°/5°/10°/15°	CS 2929-00/-05/-10/-15
obelisc end piece, oval, bronze, 32 x 26 mm		0°/5°/10°/15°/20°	CS 2930-00/-05/-10/-15/-20
Optional Implants		Angle	Prod. No.
obelisc lateral end piece, yellow, length 35 mm		0°/5°/10°/15°	CS 2914-00/-05/-10/-15
obelisc lateral end piece, green, length 45 mm		0°/5°/10°/15°	CS 2915-00/-05/-10/-15
obelisc lateral end piece, blue, length 55 mm		0°/5°/10°/15°	CS 2916-00/-05/-10/-15

Information on sterile packed implants

To order sterile packed implants, add "S" to the end of the product number (example: CS 2920-17-S).



ulrich GmbH & Co. KG
 Buchbrunnenweg 12
 89081 Ulm, Germany
 T +49 731 9654-0
 spine@ulrichmedical.com
 www.ulrichmedical.com

ulrich medical France SAS
 10, Allée Aristide Maillol
 31770 Colomiers, France
 T +33 5 34 50 91 02
 info@ulrichmedical.fr
 www.ulrichmedical.fr