

tezo™

titanium cage family



Versatile
cage family

Ulrich
medical

System

- Implant family for lumbar to lumbosacral interbody fusion
- Three cage types for posterior and anterior approaches (PLIF, TLIF, ALIF)
- Sophisticated tooting for easy insertion and high primary stability
- Large filling volume for bone/ bone substitute material



One functional principle for all cage types



Available in three sizes (S, M, L) as well as different heights and angles



Advantages

- Inserter with the same functional principle for all cages
- Rough surface structure on the inside to promote osteo-integration
- Sterile packaging – 9½ years shelf life



Components

Implants	Length x width	Height	Angle	Product number
tezo-P	24x10 mm	7 to 13 mm	5°	CS 3401-07 bis -13
tezo-P	29x10 mm	7 to 13 mm	5°	CS 3402-07 bis -13
tezo-P	29x10 mm	8 to 13 mm	12°	CS 3403-08 bis -13
tezo-T	29x11 mm	7 to 13 mm	5°	CS 3415-07 bis -13
tezo-T	29x14 mm	7 to 13 mm	5°	CS 3416-07 bis -13
tezo-T	34x16 mm	7 to 13 mm	5°	CS 3417-07 bis -13
tezo-A	28x22 mm	8 to 16 mm	8°	CS 3431-08 bis -16
tezo-A	31x25 mm	8 to 16 mm	0°	CS 3432-08 bis -16
tezo-A	31x25 mm	8 to 16 mm	8°	CS 3433-08 bis -16
tezo-A	31x25 mm	8 to 16 mm	12°	CS 3434-08 bis -16
tezo-A	34x28 mm	8 to 16 mm	0°	CS 3435-08 bis -16
tezo-A	34x28 mm	8 to 16 mm	8°	CS 3436-08 bis -16
tezo-A	34x28 mm	10 to 16 mm	12°	CS 3437-10 bis -16

Optional

Implants	Length x width	Height	Angle	Product number
tezo-P	24x10 mm	14 mm	5°	CS 3401-14
tezo-P	29x10 mm	14 mm	5°	CS 3402-14
tezo-P	29x10 mm	14 mm	12°	CS 3403-14
tezo-T	29x11 mm	14 mm	5°	CS 3415-14
tezo-T	29x14 mm	14 mm	5°	CS 3416-14
tezo-T	34x16 mm	14 mm	5°	CS 3417-14
tezo-A	28x22 mm	18 mm	8°	CS 3431-18
tezo-A	31x25 mm	18 mm	0°	CS 3432-18
tezo-A	31x25 mm	18 mm	8°	CS 3433-18
tezo-A	31x25 mm	18 mm	12°	CS 3434-18
tezo-A	34x28 mm	18 mm	0°	CS 3435-18
tezo-A	34x28 mm	18 mm	8°	CS 3436-18
tezo-A	34x28 mm	18 mm	12°	CS 3437-18

patented
or/and
pat. pend.

Over a Century
of Innovation

ulrich
medical

ulrich GmbH & Co. KG | Buchbrunnweg 12 | 89081 Ulm | Germany
Phone: +49 (0)731 9654-0 | Fax: +49 (0)731 9654-2705
spine@ulrichmedical.com | www.ulrichmedical.com