



ADDplus™

anterior distraction device plus



System

- Distractible vertebral body replacement with solid plate for the anterior cervical and cervico-thoracic spine
- Increasing angles with increasing implant size (0°–18°)
- Optionally cancellous bone screws or expandable screws
- Large filling volume for bone/ bone substitute material





Advantages

- Vertebral body substitute and plate, in a single system
- Optimized implant anchoring by expandable screw
- Smoothly distractable for individual anatomical adaptation
- Large filling volume



Components

Implants	Product number
ADDplus, Ø 12 mm, angle 0°, Height 13 – 18 mm	CS 2253-12-13
ADDplus, Ø 12 mm, angle 6°, Height 17 – 26 mm	CS 2253-12-17
ADDplus, Ø 12 mm, angle 12°, Height 25 – 41 mm	CS 2253-12-25
ADDplus, Ø 12 mm, angle 18°, Height 40 – 65 mm	CS 2253-12-40
osmium screw incl. bolt, gold, Ø 5.0 mm, Length 14 mm	CS 1300-14T
osmium screw incl. bolt, green, Ø 5.0 mm, Length 16 mm	CS 1300-16T
osmium screw incl. bolt, blue, Ø 5.0 mm, Length 18 mm	CS 1300-18T
Cancellous bone screw, Ø 4.0 mm, Length 12 – 26 mm (2 mm steps)	CS 1301-12T – CS 1301-26T
Locking screw	CS 2259

Optional

osmium screw incl. bolt, grey, Ø 4.0 mm, length 12 mm	CS 1303-12
osmium screw incl. bolt, gold, Ø 4,0 mm, length 14 mm	CS 1303-14
osmium screw incl. bolt, green, Ø 4,0 mm, length 16 mm	CS 1303-16
osmium screw incl. bolt, blue, Ø 4,0 mm, length 18 mm	CS 1303-18

Over a Century
of Innovation



ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany
Phone: +49 (0)731 9654-0 | Fax: +49 (0)731 9654-2705
spine@ulrichmedical.com | www.ulrichmedical.com